

PLATINUM GMT™ ELECTROFITNESS SYSTEM CLOUD

 **Global
Medical
Technology**
www.gmt-europa.com

C/ Isla de Alegranza 2, Nave 51
28703 San Sebastián de los Reyes
Madrid _ España
phone: (+34) 902 02 41 91
fax.: (+34) 916 53 50 19
email: Info@Gmt-Europa.com
skype: Formacion.gmt

PLATINUM GMT™ ELECTROFITNESS SYSTEM CLOUD



SPECIFICATIONS

- 8"- TOUCH COLOR SCREEN
- 4 m SPIRAL CABLE
- 11 INDEPENDENT CHANNELS (ELECTRODE PAIRS)
- 6 PRE-DEFINED TREATMENTS
- 3 CONFIGURABLE TREATMENTS
- GRADUAL CONTROL OF:
 - GENERAL INTENSITY AND EACH CHANNEL
 - CONTRACTION AND RELAXATION TIME
 - ASCENDING CONTRACTION RAMP
 - PULSE FREQUENCY
 - PULSE WIDTH
- PROGRAMMABLE REGRESSION TIME COUNTER
- PROGRESSIVE TIME COUNTER
- IDENTIFICATION CHIP CARDS FOR CLIENT AND OPERATOR
- GMT CLOUD SYSTEM
- SPECIAL TRAINING VEST:
 - NEOPRENE, ANTIBACTERIAL
 - VELCRO CLOSURE
 - ADJUSTABLE WITH FLEXIBLE STRAPS
 - 11 PAIRS OF RUBBER ELECTRODES
 - HIGH QUALITY METAL CONNECTORS

SUPPLIES

- PLATINUM GMT™ ELECTROFITNESS SYSTEM CLOUD
- NEOPRENE VEST WITH 22 ELECTRODES
- WI-FI CONNECTION ACCESSORIES FOR GMT CLOUD SYSTEM
- USER MANUAL AND WORK PROTOCOL
- TRAINING COURSE AND CERTIFICATES

APPLICATIONS

EXERCISE: Strengthens and tones muscles, preventing strain to the joints and heart, six times faster than weight training. Corrects muscle and postural imbalance. Improves physical shape with noticeable results in ten sessions.

HEALTH: Alleviates back and joint pain. Improves posture and stimulates the metabolism. Improves blood flow. Prevents osteoporosis. Releases endorphins, induces relaxation, and the muscle massage will increase well-being.


BEAUTY: Weight-loss, decrease in size, and improves cellulite. Reduces flaccidity by toning muscles and resulting in smoother tissue.

PLATINUM GMT™ ELECTROFITNESS SYSTEM CLOUD IDEAL FOR A BODY IN SHAPE

Platinum GMT™ ELECTROFITNESS System Cloud uses MES technology (muscular electrostimulation), a new and effective way to work the entire body with **muscular stimulation using mild electric impulses**. Training with this machine simultaneously stimulates the most relevant groups of superficial and core musculature in **sessions that last only 20 minutes**, which would equal a total of 10 hours of traditional exercise.

The **special neoprene antibacterial training vest** assures adequate hygiene levels, maintaining **the electrodes in the correct position**. It has velcro closures so that the electrodes can be placed in the desired position. Both the vest and the cables are equipped with **high quality metal cables**. The vest does not limit movement and the **size is adjustable** due to the flexible straps. It consists of 11 **pairs of rubber electrodes** (each one moving a muscle group), digitally controlled with the possibility of controlling contraction time, relaxation, ramp, frequency, and pulse width.

The device may be used for exercise, rehabilitation, for preventative exercise, or for esthetic purposes, where the **results can be seen in a short amount of time**. In addition, it features a remote connection through the Internet for productivity and technical monitoring (**GMT Cloud system**).



**CLOUD GMT
REMOTE
CONTROL
SYSTEM**

The **Platinum GMT™** Series features the exclusive Cloud technology with **remote control system**. Cloud GMT incorporates **operation transmission data** from the devices to the cloud (Cloud Computing) for a **complete centres management**. This allows the business owner to access the system from any Internet connected gadget (computer, tablet, Smartphone) and see in real time the work being done at the centre without the possibility of errors or data manipulation as it accesses the machine directly.

The **Cloud GMT** system also allows our **customer service** to carry out **preventive checks** that may **troubleshoot malfunctioning devices remotely** or anticipate computer updates, **reducing downtime due to repairs**. Having real-time warnings of any incident or malfunctioning can be taken care of immediately.

Cloud GMT uses a **web application** developed exclusively to manage each device. The user can access the device with the appropriate passwords and **download any information to a Internet connected gadget**.


In addition to **automatically generating and sending periodic reports** of each session; it allows the user **to setup any alarms** regarding different device performance indicators such as: power, temperature, frequency and others. It can also **generate information about the clinic** use of the devices (number of pulses per area, sessions per day, etc.), sending warnings to responsible parties if the parameters are outside set ranges. This information comes to the owner's cell phone automatically, allowing you to have **complete control of your business from anywhere in the world at any time**.

- > Maximum / minimum billings.
- > Information about workers and customers.
- > Excess time used by an operator in a session.
- > Turning on/off the device when it should not be in use.
- > Ability to prepare custom-made reports.

GMT is the only company that offers this technology.

**CLOUD GMT
REPORT SAMPLE**

| CENTRE: USERGMT | | | DEVICE: E-09999-13-13 | | | | DEVICE TYPE: ELECTRO BODY SYSTEM | | | |
|-----------------|-----------|-----------------|-----------------------|----------|----------|--------------|----------------------------------|----------------|-------------|------------|
| Nº | DATE | TREATMENT | START TIME | END TIME | DURATION | MINUTES DONE | TOTAL MINUTES | TOTAL SESSIONS | OPERATOR | CLIENT |
| 1 | 04/04/201 | CELLULITE | 9:10:24 | 9:35:59 | 25 MIN. | 25 | 565 | 29 | 00000F051FF | 000768BF2F |
| 2 | 04/04/201 | MUSCLE TONING 1 | 9:53:16 | 10:27:31 | 32 MIN. | 25 | 531 | 28 | 00000F051FF | 000768BF2E |
| 3 | 04/04/201 | MUSCLE TONING 1 | 11:21:58 | 11:43:31 | 22 MIN. | 20 | 511 | 27 | 00000F051FF | 000768C0C1 |
| 4 | 04/04/201 | MUSCLE TONING 2 | 12:15:16 | 12:47:56 | 26 MIN. | 25 | 486 | 26 | 00000F051FF | 000768BF2C |
| 5 | 03/04/201 | MUSCLE TONING 1 | 9:15:24 | 9:40:59 | 25 MIN. | 25 | 461 | 25 | 00000F051FF | 000768C151 |
| 6 | 03/04/201 | FAT-LOSS 1 | 09:53:16 | 10:23:31 | 27 MIN. | 25 | 436 | 24 | 00000F051FF | 000768BF57 |
| 7 | 03/04/201 | MUSCLE TONING 1 | 11:21:58 | 11:33:31 | 12 MIN. | 10 | 426 | 23 | 00000F051FF | 000768C0C0 |
| 8 | 03/04/201 | MUSCLE TONING 1 | 12:15:16 | 12:47:56 | 25 MIN. | 25 | 401 | 22 | 00000F051FF | 000768C11D |
| 9 | 03/04/201 | MUSCLE TONING 2 | 12:30:38 | 12:56:24 | 26 MIN. | 25 | 476 | 21 | 00000F051FF | 000768BFFD |
| 10 | 03/04/201 | CELLULITE | 15:16:37 | 15:38:12 | 26 MIN. | 25 | 451 | 20 | 00000F051FF | 00074D5B80 |



GLOBAL MEDICAL TECHNOLOGY

COMPANY

Global Medical Technology, LLC (GMT) began operations in 1998 in the field of aesthetics, specializing in the **development, marketing and maintenance of equipment intended for hair removal, weight loss and skin treatments using non-invasive methods.**

GMT manufactures and distributes **leading aesthetics devices**, according to each country's standards. GMT represents a **quality benchmark at a competitive cost** distinguished by **unique developed systems.**

The constant investment in Research and Development makes the **Platinum GMT™ Series** one of the **most well-known and technologically advanced devices** in the world. Over 1,500 cosmetic clinics in different countries use GMT devices.

GMT COMPETITIVE ADVANTAGES

- > **Proven Results.** Many years on the market and the large number of clinics using our technology are the main guarantee of quality.
- > **Exclusive GMT Cloud technology** allows owners full remote control of the treatments at each clinic (sessions, schedules, operators...), accessing the system in real time from anywhere via Internet-connected devices. It also allows GMT customer service to monitor the equipment. This helps prevent possible breakdowns and facilitates online repair when needed, reducing downtime of the device.
- > **The GMT Department of Education trains and certifies all personnel** in charge of providing aesthetic services, and answers treatments related questions. This ensures that no competitor will have the same level of service, knowledge or attention to their customers.
- > **Treatment Speed.** Compared to similar devices currently in the market, GMT's have the highest performance pulses per minute and treatments spot size, ensuring the most efficient service.
- > **High performance.** Its unique cooling system with sapphire crystals allows these devices to be working for hours without overheating.
- > **Easy to use.** Its proprietary and intuitive software automatically calculates all parameters in just three steps.
- > **Computer Monitoring.** To keep both users and clients records.
- > **System Upgrades.** Allows clinics using these devices to ensure appropriate technological responses to competitors' improvements.

TECHNICAL SERVICES

- > Ongoing training and counseling
- > Cloud GMT Service: remote control of the devices
- > Warranty
- > Sale of parts and accessories guaranteed
- > Technical Services Contracts
- > GMT-business owner coordinated operation of the devices
- > Maintenance and infrequent repairs

Global Medical Technology has its own **technical service** so that in case of breakdown customer response is immediate and in most cases the problem fixed in less than 24 hours after.

